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³ Zimmer Research Report_WA_2179_10_Rev. 1

ASTM F2052: Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Environment;

ASTM F2119: Standard Test Method for Evaluation of MR Image Artifacts from Passive Implants;

ASTM F2182: Standard Test Method for Measurement of Radio Frequency Induced Heating on or Near Passive Implants During Magnetic Resonance Imaging;

ASTM F2213: Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment;

ASTM F2503: Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment;

Catálogo de produtos MEDENS 2018 <https://www.MEDENS.com.br/arquivos/ptbr/Catalogo-de-Produtos/Catalogo-MEDENS-2018.pdf>

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WEBSTER, J. G.; 1988. Encyclopedia of Medical Devices and Instrumentation, John Wiley & Sons. v.3

SAFETY INFORMATION IN MAGNETIC RESONANCE ENVIRONMENT (RM)**INTRODUCTION:**

Emfils evaluated the behavior of dental implants and secondary components when subjected to the magnetic fields of magnetic resonance imaging (MRI) imaging studies, studies and publications through equivalence and comparison with our products. This document summarizes the results of these studies. It is intended to serve as a guide for the health professional to evaluate whether patients with Emfils implants and restorations should or should not undergo MRI.

GLOSSARY:

Gauss (G): unit of magnetic flux density. One Gauss equals 10⁻⁴ tesla;

MRI: Magnetic resonance imaging;

mN: Millinewton (N / 1000);

mT / cm: MilliTesla per centimeter (T / 1000);

RF: Radio frequency;

MRI: Magnetic resonance imaging;

SAR: SPECIFIC ABSORPTION RATE. Specific Absorption Rate (SAR) is a measure of the rate at which energy is absorbed by the human body when exposed to a radio frequency (RF) electromagnetic field. It is defined as the power absorbed by mass of fabric and has units of watts per kilogram (W / kg);

SE: Spin-echo. The spin-echo pulse sequences are one of the earliest developed and the most widely used of all magnetic resonance pulse sequences. The pulse sequence timing can be adjusted to provide T1-weighted, T2-weighted, and proton-weighted images;

T (Tesla): unit for magnetic flux density;

LITERARY REVIEW:

Some literature cited in this work is old, but the physics of the process did not change.

MRI (Magnetic Resonance Imaging) is a diagnostic tool used by healthcare professionals to produce cross-sectional images of a patient's anatomy. This technology uses magnetic fields and radio waves to create body images. Because dental implants and prostheses are made of metal or include metal parts, patients and dental professionals are often asked questions about the safety of such elements in the high-intensity magnetic fields used in the MRI environment. According to WEBSTER (1988), major concerns include magnetically induced displacement forces, radio frequency (RF) induced torque and heating effects, and affect patient safety and image artifacts that may interfere with image quality.

Emfils implants are manufactured in Grade 2 titanium (Ti) or titanium alloy (Ti6Al4V) grade 5. According to MACEDO (2008), these materials are considered paramagnetic and, therefore, their interaction with the magnetic fields is very subtle.

Patients who have metal devices, materials and implants due to paramagnetic properties can safely undergo an MRI procedure. (2001).

Emfils secondary components, such as abutments, temporary fixation screws and healing screws, are manufactured from Grade 5 titanium alloy (ASTM F136). According to MACEDO (2008), this material is considered paramagnetic or non-magnetic.

There is also the commercialization of components in a special alloy of cobalt-chromium-molybdenum (ASTM F1537), designed to the manufacture of items denominated as calcinables or UCLAS, that undergo the process of overfilling. According to SANTOS (2012) CoCrMo is a high temperature alloy frequently used in the production of prostheses and crowns, it is generally considered as non-magnetic.

The term "generally" is used because the CoCrMo alloy, in its natural state, has magnetic properties. According to SANTOS (2012), when subjected to

temperature above 1075 °C, the CoCrMo alloy loses ferromagnetism and becomes paramagnetic through allotropic transformation.

The components of CoCrMo Emfils must undergo melting, with a temperature of up to 1250 °C, that is, heating the part is inherent to its use, transforming it into paramagnetic and consequently harmless in the MRI environment.

In addition to the literature review, we also performed an equivalence analysis with the Emfils products through 3 published papers that evaluated the behavior of dental implants and components in the MRI environment meeting the standards of ASTM F 2503 - Standard for Marking Medical Devices and Other Safety Items in the MRI Environment:

¹MAX SAFETY ASSURANCE IN ENVIRONMENT OF VARIOUS MATERIALS IN TEST SAMPLES (MAGNETIC DISPLACEMENT AND ARTEFATOS)

MAGNETIC SHIFT: This study evaluated the effects of displacement forces on the main static magnetic field of a 3 Tesla MRI equipment according to ASTM F2052-14, as well as the calculation of the expected magnetic force on the magnetic field gradient more intense than 300 mT / cm, according to ASTM F2213-06 (2011) when subjected to dental materials such as Titanium (Ti), Titanium alloys (Ti6Al4V) and Cobalt chromium alloys (CoCrMo).

Therefore, the Test revealed displacement forces of less than 100 mN (in a static magnetic field of 3 T and with a mean total body SAR of 2 W / kg in normal operation).

IMAGE ARTIFACTS: Image artifacts obtained from 3 Tesla were also evaluated according to ASTM F2119-07 (2013). In a 3T RM system, the approximate length of the image artifacts for complete arcade mounts was:

- 26 mm (Ti) or 42 mm (CoCrMo) when submitted to image acquisition using gradient echo sequences;

- 13 mm (Ti) or 20 mm (CoCrMo) when submitted to image acquisition using SE sequences, respectively;

Therefore, the emergence of artifacts is expected and should be considered when necessary the analysis of images. Image artifacts pose no risk to the patient. In the analysis of the image the professional will verify and circumvent the possibly distorted images.

² ASSURANCE OF SAFETY IN ENVIRONMENT OF RM OF CONSTRUCTIONS WITH DENTAL IMPLANTS (MAGNETIC SHIFTING, HEATING BY RF AND ARTEFATOS)

RF HEATING: The evaluation of the effects of heating caused by the radiofrequency field (RF field) during MR exams with 1.5 Tesla and 3 Tesla fields according to ASTM F2182-11a was performed. The study revealed that metallic components made of titanium (Ti), titanium alloys (Ti6Al4V) and chromium alloys (CoCrMo) are expected to produce a maximum temperature increase of less than 5.9 °C when subjected to fields of 3 T and less than 5.2 °C when subjected to 1.5 T fields (after 15 minutes under continuous examination).

Therefore, RF heating produces a maximum temperature rise of less than 6 °C, causing low impact when considering the cooling effects generated in the adjacent tissues or by the blood circulation.

³MAX SAFETY ASSURANCE IN MRI ENVIRONMENT WITH DENTAL IMPLANTS (MAGNETIC SHIFT)

Tests have shown that while each metal material exhibits a small but measurable magnetic attraction in the Tesla 3.0 and Tesla environments of 1.5, the maximum magnetic force exerted on a device is less than 25% of the force exerted on the device due to gravity. None of the metal materials exhibited any torque movement in 1.5 and 3.0 Tesla environments.

Therefore, no movement or deflection of devices fabricated from the above-mentioned metallic materials is expected in 1.5 and 3.0 Tesla environments.

CONCLUSION:

Emfils products have material, volume, functional and geometric

compatibility with the products tested in the works cited above. The results of these studies indicate that MEDENS original implants and components manufactured from Titanium (Ti), Titanium alloys (Ti6Al4V) and CoCrMo alloys can be safely used in the MRI environment considered as “conditional MRI” according to standard ASTM F2503.

“Conditional MRI” indicates that Emfils implants and components are safe to the user as long as they meet the following conditions:

- Static magnetic field of 3 and 1.5 Tesla only.
- Spatial gradient field of 720 Gauss / cm, or less (0.072 Tesla / centimeter).
- The maximum magnetic resonance system reported average spectral absorption rate (SAR) of 2W / kg per 15 minutes of scanning.

Under the conditions of use recited above, the product is expected to produce a temperature below 6 ° after 15 minutes of continuous exposure, being considered low impact and harmless.

The image artifact caused extends to 26mm for titanium components and 42mm for CoCrMo alloy when exposed to a sequence of echo-gradient pulses in a 3 Tesla magnetic resonance environment, and does not compromise the health of the patient.

DISCLAIMER OF LIABILITY:

Test results indicate that Emfils Implants manufactured from the materials listed above can be considered safe for use in the “conditional” MRI environment, ie in the testing environment. However, Emfils assumes no responsibility for its use in the clinical RM environment. Emfils recommends that medical professionals and the manufacturer of MRI equipment be consulted for compatibility with implants before use.

Products considered “Conditional MR” are identified by the MR symbol on the label.

If there is no MR symbol on the product label, please note that the product has not been evaluated for safety and compatibility in the MR environment. The product has not been tested for heating or displacement in the MR environment.

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³Zimmer Research Report_WA_2179_10_Rev. 1

ASTM F2052: Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Environment;

ASTM F2119: Standard Test Method for Evaluation of MR Image Artifacts from Passive Implants;

ASTM F2182: Standard Test Method for Measurement of Radio Frequency Induced Heating On or Near Passive Implants During Magnetic Resonance Imaging;

ASTM F2213: Standard Test Method for Measurement of Magnetically Induced Torque on Medical Devices in the Magnetic Resonance Environment;

ASTM F2503: Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment;

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